

Usage instructions:

PLEASE READ PRIOR TO USE

CHARGING

- Fully charge battery for up to 4 hours before first use from the device.
- To charge, connect the DC end of the supplied 12.6V adapter to the battery's charging port and plug AC side cable into wall outlet.
- On adaptor red lights pulsate to show indicate active charging. On adaptor green lights pulsate to show is fully charged.
- With battery attached. LED screen on unit to show battery level.
- Battery may be recharged at any time and at any battery level.
- It is not recommended to fully deplete the battery below 20% level.
- The average run time is 3 hours depending on speed level and pressure applied during use.

SAFETY INSTRUCTIONS

DANGER: TO REDUCE RISKS OF ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY; THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

- For adult use only.
- Use only on a dry clean surface of the body over clothing by lightly pressing and moving across the body for approximately 60 seconds per region.
- Use the massage gun only on the body's soft tissue as desired without producing pain or discomfort. Do not use on head or any hard or bony area of the body.
- Use only the applicators that provide the best desired effect.
- Bruising can occur regardless of control setting or pressure applied. Check treated areas frequently and immediately stop at first sign of pain or discomfort.
- Keep fingers, hair or any other body part away from the shaft and backside of the applicator as pinching may occur.
- Do not place any objects into the ventilation holes of the Massage Gun
- Do not immerse in water or allow water to enter the ventilation holes of the Massage Gun
- Do not drop or misuse the Massage Gun.
- Only recharge with supplied 12.6VDC Massage Gun charger.
- Carefully examine device and battery before each use.
- Do not tamper or alter Massage Gun in any way.
- Never leave the Massage Gun operating or charging unattended.

Disclaimer

WARNINGS & CAUTION

- Do not immerse in water. Keep away from liquids or heat sources. Keep ventilation ports free from dust and debris.
- Do not remove screws or attempt to disassemble.
- Do not operate continuously for more than one hour. Allow device to rest for 30 minutes before re-using.
- Unplug the charger from device battery after charging or prior to use.
- For adult use only. Do not use if injured. Consult your doctor before using this product.

SAFETY INSTRUCTIONS

DANGER: TO REDUCE RISKS OF ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY; THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

- For adult use only.
- Use only on a dry clean surface of the body over clothing by lightly pressing and moving across the body for approximately 60 seconds per region.
- Use the massage gun only on the body's soft tissue as desired without producing pain or discomfort. Do not use on head or any hard or bony area of the body.
- Use only the applicators that provide the best desired effect.
- Bruising can occur regardless of control setting or pressure applied. Check treated areas frequently and immediately stop at first sign of pain or discomfort.
- Keep fingers, hair or any other body part away from the shaft and backside of the applicator as pinching may occur.
- Do not place any objects into the ventilation holes of the Massage Gun
- Do not immerse in water or allow water to enter the ventilation holes of the Massage Gun
- Do not drop or misuse the Massage Gun.
- Only recharge with supplied 12.6VDC Massage Gun charger.
- Carefully examine device and battery before each use.
- Do not tamper or alter Massage Gun in any way.
- Never leave the Massage Gun operating or charging unattended.

Disclaimer

PLEASE DO NOT USE THE MESSAGE GUN OR ANY VIBRATION DEVICE WITHOUT FIRST OBTAINING APPROVAL FROM YOUR DOCTOR IF ANY OF THE FOLLOWING APPLY:

Pregnancy, diabetes with complications such as neuropathy or retinal damage, wear of pace-makers, recent surgery, epilepsy or migraines, herniated disks, spondylolisthesis,

spondylo-ly-sis, or spondylosis, recent joint replacements or IUD's, metal pins or plates or any concerns about your physical health. Frail individuals and children should be accompanied by an adult when using any vibration device. These contra indications do not mean that you are not able to use a vibration or exercise device but we advise you to consult a doctor first.

Ongoing research is being done to the effects of vibration exercise for specific medical disorders. This will likely result in the shortening of the list of contra indications as shown above. Practical experience has shown that integration of vibration exercise into a treatment plan is advisory in a number of cases. This must be done on the advice of and in accompaniment with a doctor, specialist or physiotherapist.

WARRANTY

IF WITHIN ONE YEAR FROM THE DATE OF PURCHASE, THIS PRODUCT FAILS DUE TO A DEFECT IN MATERIAL OR WORKMANSHIP WIU REPAIR OR REPLACE THE PRODUCT, OR NECESSARY COMPONENTS, FREE OF CHARGE.

THIS WARRANTY EXCLUDES:

- (a) Damage caused by accident, abuse, mishandling, or transport.
- (b) Units subjected to unauthorized repair.
- (c) Units not used in accordance with instructions.
- (d) Damage exceeding the cost of the product.
- (e) Deterioration of the delivered product resulting from abnormal storage and/or safeguarding conditions on the client's premises.
- (f) Failure to provide the dated proof of purchase.

Some states and countries do not allow a limitation of damages, so the foregoing limitation may not apply to you. This warranty guarantees specific legal rights, and but other rights may vary from country to country and from state to state. This warranty service is available through the operating policies and international distributors.

Simply Brands

www.simply-brands.com




VIBRATION MASSAGE DEVICE

Instruction Manual



Packing list:



Model: 20w

Working Time: 1-3H

N.W.: 850g

Power Supply: Non-removable lithium battery

Charging Time: 3H

Input Voltage: 100-240v

Battery Voltage: 12.6v

Importer Name:
Simply Brands (Asia) Ltd

Model Number:
SBA-00099

Importer Address:
Unit 04-05, 16th Floor, The Broadway No.
54-62 Lockhart Road, Wanchai, Hong
Kong

Please retain this information for
future reference.

2

Usage instructions:

PLEASE READ TO USE



Battery Power Display
Turn On & Gear Display
Working

Over shift “ + ”

Downshift “ - ”

Power Switch

Battery Display Lights

Gear Display Lights


Battery Power Display
Lights

Battery Display Lights

3

Usage instructions:

PLEASE READ TO USE




Charging
port

ON/OFF
Switch

4


Usage instructions:

PLEASE READ TO USE



Round Head

For large muscle groups
like quads, glutes, extreme
muscles, and joint soreness.




U-Shaped Head

For spine, neck, and shoulder



Bullet Head

For joints, deep tissues,
trigger points feet, and wrist.



Fat Head

For all parts of the body.

5